


ORTHOSCREW



- Self-perforating and self-threading micro Implants
- Gr 5 medical titanium
- Not treated surface
- Crosswise cavity: .022" x .028"
- Hex head: 3 x h 1,4 mm

- Temporary disposable micro implants.
- Minimal risk of osseointegration.
- Hexagonal, smooth rounded head for an optimum comfort.
- They allow particularly complex movements, impossible with a conventional approach.

Ø	Length	Neck Height 0,5	Neck Height 1,5
		Code	Code
1,6	5,0	0110S1605/3	01101605/3
	7,0	0110S1607/3	01101607/3
	9,0	0110S1609/3	01101609/3
	11,0	0110S1611/3	01101611/3

Packaging

- Packaging in compliance with ISO 11607-1 and 2.
- Sterilization by gamma rays 25 kGy
- Sterility guaranteed for 5 years, thanks to the double packaging
- The packaging contains:
3 implants



Ø 1,6

Body: Ø 1,65 mm

Hole: Ø 0,8 mm

Thread pitch: 0,95 mm

Tapered apex



Neck
h 0,5 mm



Length	Code
5,0	0110S1605/3
7,0	0110S1607/3
9,0	0110S1609/3
11,0	0110S1611/3

Neck
h 1,5 mm



Length	Code
5,0	01101605/3
7,0	01101607/3
9,0	01101609/3
11,0	01101611/3

Surgical instruments for ORTHOSCREW implants

Drills



Steel
Ø 1,0 mm

01DS10S



Steel
Ø 1,2 mm

11DS12S

Connectors



Manual

01TWN25



Contra-angle

01TLCAN25

Accessories



Manual
digma

01TLM



Screwdriver

01TLO



Empty box

01V71